

# Abbott Laboratories recalls blood glucose meters

April 15 2013

---

Abbott Laboratories is recalling its FreeStyle InsuLinx Blood Glucose Meters after finding that they display and store incorrect test results for dangerously high blood sugar levels.

The Abbott Park, Ill., company says the meters will display and store readings for [blood glucose levels](#) of 1,024 milligrams per deciliter or higher at much lower levels. That's an extremely high and rare level that requires immediate medical attention.

The inaccurate readings can lead to a delay in identifying the problem or incorrect treatment. Abbott says that can lead to serious injury or death.

Abbott says customers with the meter should contact its diabetes care customer service or download a software update that resolves the issue.

Johnson & Johnson announced a similar recall last month for several types of its blood glucose meters.

Copyright 2013 The Associated Press. All rights reserved. This material may not be published, broadcast, rewritten or redistributed.

Citation: Abbott Laboratories recalls blood glucose meters (2013, April 15) retrieved 3 May 2024 from <https://medicalxpress.com/news/2013-04-abbott-laboratories-recalls-blood-glucose.html>

This document is subject to copyright. Apart from any fair dealing for the purpose of private
---

study or research, no part may be reproduced without the written permission. The content is provided for information purposes only.