


# Software visualizes the heart for doctors

December 9 2013

---



Siemens developed a procedure for angioplasty which supports doctors to make interventions at the heart faster and more precise. With angioplasty doctors open up an artery with a stent, potentially staving off more invasive procedures like heart bypass surgery. To deploy the stent, doctors follow its position by X-ray fluoroscopy, but the procedure is limited by poor visibility and difficult use. This can be helped by use of intravascular ultrasound (IVUS) which implies the use of costly catheters. The visualizing tool CLEARstent makes this procedure easier. The picture shows a stent during the procedure.

Siemens developed a procedure for angioplasty which supports doctors to make interventions at the heart faster and more precise. With angioplasty doctors open up an artery, potentially staving off more invasive procedures like heart bypass surgery. This method positions and

expands a balloon, enclosed in a stent a metal stent, to open up a region in a coronary artery that can clog. To deploy the stent, doctors follow its position by X-ray fluoroscopy, but the procedure is limited by poor visibility and difficult use. This can be helped by use of intravascular ultrasound (IVUS) which implies the use of costly catheters.

The visualizing tool CLEARstent makes this procedure easier. The software is for several angiography devices from Siemens which all work with low doses of radiation. After a few seconds imaging acquisition time, or using previously acquired images, CLEARstent can tell clinicians if the stent is properly positioned and expanded. Doctors need to figure out where they want to place the stent and afterwards whether the stent is fully embedded in the blood vessel's wall. CLEARstent enables both these processes.

Experts from Siemens Healthcare and the global Siemens research organization Corporate Technology (CT) in Princeton, New Jersey, developed the software. CT researchers developed also the software CLEARstent Live which improves visibility, providing better images during live angioplasty when the heart is beating. Imaging of the stent by traditional X-ray fluoroscopy is complicated by the motion induced by the pumping heart. X-ray fluoroscopy tends to pick up blurry, faint images of the stent.



X-ray fluoroscopy tends to pick up blurry, faint images of the stent.

To improve the visibility of the stent, the CT team developed a robust method that detects the balloon markers at the stent's ends by a learning step with early snapshots of the stent. This preparative step would help the program identify the stent during live image acquisitions. To reduce blurry images further, CLEARstent Live performs a common method known as image registration. This method deforms and aligns the clearly defined [stents](#) in adjacent frames based on the marker position, thereby stabilizing the stent during real-time imaging. Another benefit for the doctor is that CLEARstent Live can overlay the [images](#) of the stent with those of the surrounding blood vessels. This ensures that live [angioplasty](#) is performed under the most physiological condition.

Provided by Siemens

Citation: Software visualizes the heart for doctors (2013, December 9) retrieved 6 May 2024 from <https://medicalxpress.com/news/2013-12-software-visualizes-heart-doctors.html>

This document is subject to copyright. Apart from any fair dealing for the purpose of private study or research, no part may be reproduced without the written permission. The content is provided for information purposes only.