

May be room for improvement in U/S transducer hygiene

January 23 2015



(HealthDay)—For endoluminal procedures relying on barrier protection to avoid contamination, permeability of materials may not always be considered, according to a study published in the February issue of the *Journal of Clinical Ultrasound*.

Mikael Häggström, M.D., from Uppsala University in Sweden, and colleagues examined whether current <u>hygiene practices</u> are appropriate during sonographic examinations. Questionnaires were completed for 25 departments in five major hospitals in Sweden by personnel who were responsible for or acquainted with the local hygiene procedures.

The researchers found that for transvaginal and transrectal sonographic examinations, the most common method for decontamination of the transducer was barrier protection during the procedure and alcohol



cleansing after the procedure. The predominant cover material was latex. One department used polyethylene gloves and another used nitrile gloves—both for transvaginal ultrasonography. All hospitals were using alcohol and paper or cloth for decontamination at a minimum in transcutaneous examinations. Barrier protection was not used in transesophageal examinations, and <u>decontamination</u> was performed with an alkylating agent.

"The hygiene practices appear to be appropriate at most hospitals, but there is a prevalence of transducer cover materials of unacceptable <u>permeability</u>, as well as use of gloves on transducers despite insufficient evidence of safety," the authors write.

More information: Abstract

Full Text (subscription or payment may be required)

Copyright © 2015 HealthDay. All rights reserved.

Citation: May be room for improvement in U/S transducer hygiene (2015, January 23) retrieved 12 May 2024 from <u>https://medicalxpress.com/news/2015-01-room-transducer-hygiene.html</u>

This document is subject to copyright. Apart from any fair dealing for the purpose of private study or research, no part may be reproduced without the written permission. The content is provided for information purposes only.