


How to identify drugs that work best for each cancer patient

April 22 2015



MIT chemical engineers have designed an implantable device that can deliver many drugs at once, allowing researchers to determine which drugs are the most effective against a patient's tumor.

More than 100 drugs have been approved to treat cancer, but predicting

which ones will help a particular patient is an inexact science at best.

A new [device](#) developed at MIT may change that. The implantable device, about the size of the grain of rice, can carry small doses of up to 30 different drugs. After implanting it in a tumor and letting the drugs diffuse into the tissue, researchers can measure how effectively each one kills the patient's cancer cells.

Such a device could eliminate much of the guesswork now involved in choosing cancer treatments, says Oliver Jonas, a postdoc at MIT's Koch Institute for Integrative Cancer Research and lead author of a paper describing the device in the April 22 issue of *Science Translational Medicine*.

"You can use it to test a patient for a range of available drugs, and pick the one that works best," Jonas says.

The paper's senior authors are Robert Langer, the David H. Koch Professor at MIT and a member of the Koch Institute, the Institute for Medical Engineering and Science, and the Department of Chemical Engineering; and Michael Cima, the David H. Koch Professor of Engineering at MIT and a member of the Koch Institute and the Department of Materials Science and Engineering.

Putting the lab in the patient

Most of the commonly used [cancer drugs](#) work by damaging DNA or otherwise interfering with cell function. Recently, scientists have also developed more targeted drugs designed to kill tumor cells that carry a specific genetic mutation. However, it is usually difficult to predict whether a particular [drug](#) will be effective in an individual patient.

In some cases, doctors extract [tumor cells](#), grow them in a lab dish, and

treat them with different drugs to see which ones are most effective. However, this process removes the cells from their natural environment, which can play an important role in how a tumor responds to drug treatment, Jonas says.

"The approach that we thought would be good to try is to essentially put the lab into the patient," he says. "It's safe and you can do all of your sensitivity testing in the native microenvironment."

The device, made from a stiff, crystalline polymer, can be implanted in a patient's tumor using a biopsy needle. After implantation, drugs seep 200 to 300 microns into the tumor, but do not overlap with each other. Any type of drug can go into the reservoir, and the researchers can formulate the drugs so that the doses that reach the [cancer cells](#) are similar to what they would receive if the drug were given by typical delivery methods such as intravenous injection.

After one day of drug exposure, the implant is removed, along with a small sample of the tumor tissue surrounding it, and the researchers analyze the drug effects by slicing up the tissue sample and staining it with antibodies that can detect markers of cell death or proliferation.

Ranking cancer drugs

To test the device, the researchers implanted it in mice that had been grafted with human prostate, breast, and melanoma tumors. These tumors are known to have varying sensitivity to different cancer drugs, and the MIT team's results corresponded to those previously seen differences.

The researchers then tested the device with a type of [breast cancer](#) known as triple negative, which lacks the three most common breast cancer markers: estrogen receptor, progesterone receptor, and Her2. This

form of cancer is particularly aggressive, and none of the drugs used against it are targeted to a specific genetic marker.

Using the device, the researchers found that triple negative tumors responded differently to five of the drugs commonly used to treat them. The most effective was paclitaxel, followed by doxorubicin, cisplatin, gemcitabine, and lapatinib. They found the same results when delivering these drugs by intravenous injection, suggesting that the device is an accurate predictor of drug sensitivity.

In this study, the researchers compared single drugs to each other, but the device could also be used to test different drug combinations by putting two or three drugs into the same reservoir, Jonas says.

"This device could help us identify the best chemotherapy agents and combinations for every tumor prior to starting systemic administration of chemotherapy, as opposed to making choices based on population-based statistics. This has been a longstanding pursuit of the oncology community and an important step toward our goal of developing precision-based cancer therapy," says Jose Baselga, chief medical officer at Memorial Sloan Kettering Cancer Center and an author of the paper.

The researchers are now working on ways to make the device easier to read while it is still inside the patient, allowing them to get results faster. They are also planning to launch a clinical trial in breast cancer patients next year.

Another possible application for this device is to guide the development and testing of new [cancer](#) drugs. Researchers could create several different variants of a promising compound and test them all at once in a small trial of human patients, allowing them to choose the best one to carry on to a larger clinical trial.

More information: "An implantable microdevice to perform high-throughput in vivo drug sensitivity testing in tumors," *Science Translational Medicine*, [stm.sciencemag.org/lookup/doi/ ... scitranslmed.3010564](https://stm.sciencemag.org/lookup/doi/10.1126/scitranslmed.3010564)

Provided by Massachusetts Institute of Technology

Citation: How to identify drugs that work best for each cancer patient (2015, April 22) retrieved 10 April 2024 from <https://medicalxpress.com/news/2015-04-drugs-cancer-patient.html>

| |
|--|
| <p>This document is subject to copyright. Apart from any fair dealing for the purpose of private study or research, no part may be reproduced without the written permission. The content is provided for information purposes only.</p> |
|--|