

New system can assess severity of von Willebrand disease

February 5 2016



(HealthDay)—A new microchip-flow chamber system (T-TAS) can be



used to discriminate and predict bleeding score (BS) in type 1 von Willebrand disease (VWD), according to a study published online Jan. 27 in the *Journal of Thrombosis and Haemostasis*.

Keiji Nogami, M.D., Ph.D., from the Nara Medical University in Japan, and colleagues examined the potential of T-TAS for assessing physiological hemostasis in VWD. Fifty type 1 VWD patients with whole blood samples were assessed with microchips coated with collagen (PL-chip) or collagen/thromboplastin (AR-chip) to examine platelet thrombus formation (PTF) at high-shear rates or fibrin-rich PTF at lowshear rates. Flow pressure curves were used to calculate the times to increase 10 and 30 kPa (T_{10} and T_{30}) from baseline. A standardized questionnaire was used to determine BS.

The researchers found that irrespective of the flow rate, PL-T₁₀ values correlated with BS better than von Willebrand factor-ristocetin cofactor activity, while there was a weak correlation for AR-T₁₀ with BS. Lower levels of VWF and higher BS were seen for patients with PL-T₁₀ >10 minutes or AR-T₁₀ >30 minutes, compared to those with PL-T₁₀ ≤10 minutes or AR-T₁₀ ≤30 minutes, with the greater difference observed with PL-T₁₀. The best correlation for clinical severity was seen with PL-T₁₀ >8 minutes.

"T-TAS could be a useful technique for discriminating and predicting the BS in VWD type 1 patients," the authors write.

One author is employed by the Fujimori Kogyo Co., the manufacturer of T-TAS.

More information: Abstract

Full Text (subscription or payment may be required)



Copyright © 2016 <u>HealthDay</u>. All rights reserved.



Citation: New system can assess severity of von Willebrand disease (2016, February 5) retrieved 27 April 2024 from

https://medicalxpress.com/news/2016-02-severity-von-willebrand-disease.html

This document is subject to copyright. Apart from any fair dealing for the purpose of private study or research, no part may be reproduced without the written permission. The content is provided for information purposes only.