

# Helmintex assay detects low intensity schistosomiasis infections

March 8 2018

---



Dusk in Candéal, Sergipe, northeastern Brazil. Credit: Graeff-Teixeira, et al,

2018

Individuals who release low numbers of schistosome eggs in their feces may not be detected by classical methods to detect schistosomiasis. Instead, the Helmintex method— which uses a magnetic field to isolate eggs—is more sensitive for these low intensity schistosomiasis infections, researchers report in *PLOS Neglected Tropical Diseases*.

Schistosomiasis is a common infection that affects over 290 million individuals, especially in Sub-Saharan Africa, Asia and South America. Classical diagnostic methods, including the Kato-Katz (KK) fecal smear recommended by the World Health Organization, are not accurate in situations where low egg burden is present in stool. Alternative methods, including a point-of-care immunodiagnostic for schistosome cathodic circulating antigen in urine (POC-CCA) and the Helmintex (HTX) [method](#), have been proposed as substitutes.

In the new work, Carlos Graeff-Teixeira of the Pontificia Universidade Catolica de Rio Grande de Sul, and colleagues compared the KK, POC-CCA and HTX methods in a population-based study. 461 children and adults from Candeal, Brazil provided stool, urine, and blood samples for the comparison; an additional 220 submitted at least one specimen type.

The prevalence of schistosomiasis was found to be 71% with the POC-CCA method, 40% with HTX, and 11% with KK. More than 70% of individuals had very low egg burdens in their feces. However, the POC-CCA method provided the lowest proportion of true positive results. The HTX method was 100% sensitive with egg burdens higher than 1.3 epg and detected [eggs](#) in 29% of the samples that were negative according to the KK method.

"The HTX method is very sensitive and should be used as a reference method for diagnosing intestinal schistosomiasis and for comparative evaluation of other tests," the researchers say. "The HTX method should also be considered for use in the monitoring and certification of transmission interruption."

**More information:** Catieli Gobetti Lindholz et al, Study of diagnostic accuracy of Helmintex, Kato-Katz, and POC-CCA methods for diagnosing intestinal schistosomiasis in Candéal, a low intensity transmission area in northeastern Brazil, *PLOS Neglected Tropical Diseases* (2018). [DOI: 10.1371/journal.pntd.0006274](https://doi.org/10.1371/journal.pntd.0006274)

Provided by Public Library of Science

Citation: Helmintex assay detects low intensity schistosomiasis infections (2018, March 8) retrieved 20 April 2024 from <https://medicalxpress.com/news/2018-03-helmintex-assay-intensity-schistosomiasis-infections.html>

This document is subject to copyright. Apart from any fair dealing for the purpose of private study or research, no part may be reproduced without the written permission. The content is provided for information purposes only.