


Stem-cell technology aids 3-D printed cartilage repair

April 20 2018, by Lea Kivivali



Fluorescent microscopy image of biospheres with stem cells growing on their surface before detachment by infra-red beam. Scale Bar = 0.1 mm. Credit: Yashaswini Vegi

Novel stem-cell technology developed at Swinburne will be used to grow the massive number of stem cells required for a new hand-held 3-D

printer that will enable surgeons to create patient-specific bone and cartilage.

The technology, called BioSphere, is being developed in collaboration with St Vincent's Hospital in Melbourne to support the development of Biopen.


The Biopen is a hand-held 3-D printer pen filled with stem cell ink that will allow surgeons to repair damaged bone and cartilage by 'drawing' new [cells](#) directly onto bone during surgery and then filling in any damaged areas.

The BioSphere technology is part of Yashaswini Vegi's Ph.D. project within the ARC Training Centre for Biodevices at Swinburne in partnership with St Vincent's.

"In the near future it will be possible to use a patient's own [stem cells](#) to quickly and successfully repair damaged cartilage and bone," says Ms Vegi's supervisor, research scientist Dr. Nicholas Reynolds.

"This is possible thanks to stem cells' unique ability to transform into different cells and tissues, enabling us to regrow bones and cartilage inside the patient.

"In order to do this, stem cells are extracted from the patient, grown in the lab and transformed into the desired tissue type before re-implantation into the injured body part via the Biopen."



Left: Biosphere coated with light sensitive gold nanoparticles that when activated cause a gentle detachment of cells grown on the spheres. Right: Biosphere with stem cells growing on its surface before detachment by infra-red beam. Biospheres are approximately 0.3 mm in diameter. Credit: Yashaswini Vegi

Dr. Reynolds says progress in this field has been slowed by a lack of means to generate the very large number of cells required for successful treatments and also a lack of methods to recover the cells from the 'dishes' they are grown on without damaging them.

The BioSphere technology will replace dishes with thousands of polymer beads less than a millimetre in diameter.

"This will vastly increase the available area for the cells to grow on enabling researchers to generate many more cells than is currently possible," Dr. Reynolds says.

"The surface of the beads will be coated with light sensitive nanoparticles that do not affect the cells' growth until the particles are

activated with an infra-red beam. At this point the activated nanoparticles cause the cells to gently detach from the beads, and can be easily recovered for re-implantation into the injured patient."

BioMedTech funding

The Biopen is one of 11 projects selected for the first round of the \$35 million BioMedTech Horizons program investment recently announced by the Minister for Health, the Hon Greg Hunt MP.

"This funding provides valuable support for the ongoing BioSphere project undertaken between Swinburne and St Vincent's Hospital Melbourne and further strengthens Swinburne's involvement in the Aikenhead Centre of Medical Discovery initiative," says Ms Vegi's co-supervisor and Bioengineering Program Leader – Iverson Health Innovation Research Institute, Professor Simon Moulton.

The BioMedTech Horizons program aims to help Australia move more cutting edge-ideas and breakthrough discoveries towards proof-of-concept and commercialisation, and stimulate collaboration across disciplines between the research, industry and technology sectors to maximize entrepreneurship and idea potential.

Provided by Swinburne University of Technology

Citation: Stem-cell technology aids 3-D printed cartilage repair (2018, April 20) retrieved 2 May 2024 from <https://medicalxpress.com/news/2018-04-stem-cell-technology-aids-d-cartilage.html>

<p>This document is subject to copyright. Apart from any fair dealing for the purpose of private study or research, no part may be reproduced without the written permission. The content is provided for information purposes only.</p>
--