


## Mapping functional connectivity in 3-D artificial brain model by analyzing neural signals

## March 10 2021



A Schematic illustrations showing three 2D multifunctional MEAs before stacking and bonding (left), assembled 3D high-density multifunctional MEA with a PDMS fluidic interface and the multifunctional shank for optical and chemical stimulations (middle), and the application to 3D neural network model in vitro compartmentalized with two somatic and neurite regions (right). Credit: Korea Institute of Science and Technology(KIST)

The human brain is less accessible than other organs because it is covered by a thick, hard skull. As a result, researchers have been limited to low-resolution imaging or analysis of brain signals measured outside the skull. This has proved to be a major hindrance in brain research,




including research on developmental stages, causes of diseases, and their treatments. Recently, studies have been performed using primary neurons from rats or human-derived induced pluripotent stem cells (iPSCs) to create artificial brain models that have been applied to investigate brain developmental processes and the causes of brain diseases. These studies are expected to play a key role to unlocking the mysteries of the brain.

In the past, artificial <u>brain</u> models were created and studied in 2D; however, in 2017, a research team from KIST developed a 3D artificial brain model that more closely resembled the real brain. Unfortunately, due to the absence of an analytical framework for studying signals in a 3D brain model, studies were limited to analyses of surface signals or had to reform the 3D structure to a flat shape. As such, tracking <u>neural signals</u> in a complex, interconnected artificial network remained a challenge.

The Korea Institute of Science and Technology (KIST) announced that the research teams of Doctors Il-Joo Cho and Nakwon Choi have developed a <u>analysis system</u> that can apply precise non-destructive stimuli to a 3D artificial neural circuit and measure neural signals in real-time from multiple locations inside the model at the cellular level.

The 3D multifunctional system for measuring neural signals is in the form of a 50µm-wide needel shaped silicon probe array (about half the width of a human hair) integrated with 63 microelectrodes. When this system is inserted in the artificial brain model, it is capable of simultaneously measuring signals from multiple locations inside the neural circuit. The probe contains an optical fiber and drug-delivery channels, enabling precise stimulation of neurons using light or drugs. By measuring functional changes in the brain model in response to these stimuli, the model can be used to study brain function and brain diseases.





A researcher of KIST is looking at a three-dimensional multifunctional electrode chip developed by Dr. Il-Joo Cho and Dr. Nakwon Choi Credit: Korea Institute of Science and Technology(KIST)

Using this system to stimulate neural circuits in the artificial brain model optically and simultaneously measure the spread of the response signal in multiple locations, the research team demonstrated that the propagation speed of neural signals were different according to directions inside the 3D brain modeln. In addition to structural brain maps, which can be constructed using electron microscopy, this study demonstrated the possibility of constructing 3D functional brain maps that show how different circuits are functionally connected within complex artificial brain networks.



Dr. Choi, from KIST, stated that, "The newly developed system allows us to study various developmental brain disorders and the causes of and treatments for brain diseases." Co-PI Dr. Cho added, "This system enables functional measurements from 3D artificial brain models, which was previously impossible. We expect that the development of this system will help to radically reduce the time required to develop drug or treatments for various brain diseases."

**More information:** Hyogeun Shin et al, 3D high-density microelectrode array with optical stimulation and drug delivery for investigating neural circuit dynamics, *Nature Communications* (2021). DOI: 10.1038/s41467-020-20763-3

## Provided by National Research Council of Science & Technology

Citation: Mapping functional connectivity in 3-D artificial brain model by analyzing neural signals (2021, March 10) retrieved 5 May 2024 from <a href="https://medicalxpress.com/news/2021-03-functional-d-artificial-brain-neural.html">https://medicalxpress.com/news/2021-03-functional-d-artificial-brain-neural.html</a>

This document is subject to copyright. Apart from any fair dealing for the purpose of private study or research, no part may be reproduced without the written permission. The content is provided for information purposes only.