


Head-mounted microscope allows long-term brain imaging in freely moving mice

December 9 2021



Credit: CC0 Public Domain

Researchers have developed the first detachable head-mounted photoacoustic microscope for imaging brain activity in freely moving mice. Because the device can be removed after imaging, it enables longterm studies that could reveal important new insights into neurodegenerative diseases and other neurological disorders.

"Epilepsy, Alzheimer's disease and Parkinson's disease can all seriously interrupt neurovascular coupling—the link between neural activity and subsequent changes in cerebral blood flow," said research team leader Lei Xi from the Southern University of Science and Technology in China. "Our new probe is ideal for studying neurovascular coupling



because it has the potential to capture the dynamics of both neuron and vascular networks simultaneously."

The new microscope probe weighs just 1.8 grams and is based on optical resolution photoacoustic microscopy (ORPAM), which can capture anatomical and functional dynamics of the brain without requiring the use of fluorescent tags or labels. In the Optica Publishing Group journal *Optics Letters*, Xi and colleagues describe how they optimized the design of the new probe to make it light enough for use in freely moving <u>mice</u>.

"Head-mounted microscopes that use multi-photon or fluorescence imaging primarily capture the activities of single neurons," said Xi. "Our ORPAM probe can capture cerebral vascular <u>network</u> and hemodynamics of large portions of the cerebral cortex with capillarylevel resolution without requiring any labels."

Miniaturizing a microscope

The new work builds on a wearable ORPAM probe the researchers previously built for freely moving rats. Although it performed well, it had to be permanently fixed on the rat and was too large and heavy to be carried around by mice, which are the preferred animal models for many brain studies.

To make a smaller probe, the researchers used optical simulation calculations to carefully optimize the entire light path of the <u>microscope</u>. They also selected high-performance miniature components including an aspheric lens, micro-electro-mechanical-system (MEMS) scanner and customized miniaturized piezoelectric ultrasonic detector.

The resulting probe weighs less than 10% of an adult mouse's body weight, has a resolution of 2.8 microns and can image a large field of view of $3\times3 \text{ mm}^2$. To make it detachable, the researchers incorporated



three pairs of magnets that connect the imaging probe with a lightweight mounting base attached to the mouse's skull. The probe can be removed directly after imaging and then reinstalled later, enabling repeated, longterm monitoring of freely moving animals.

Long-term imaging

After evaluating the probe's performance on a synthetic material that mimics soft tissue, the researchers used it to image vascular networks in the cerebral cortex of a mouse for 40 minutes. They also conducted long-term monitoring experiments that lasted for seven days. The results from these tests showed that the probe can provide stable and high-quality ORPAM <u>images</u> in freely moving mice.

"In the future, we intend to develop a probe with capillary-level resolution, video-rate imaging speed and a field of view large enough to capture the entire mouse cerebral cortex," said Xi. "We want to push the performance of the probe so that even more can be learned about brain activity and how it might relate to disease and health."

More information: Guo Heng et al, Detachable head-mounted photoacoustic microscope in freely moving mice, *Optics Letters* (2021). DOI: 10.1364/OL.444226

Provided by The Optical Society

Citation: Head-mounted microscope allows long-term brain imaging in freely moving mice (2021, December 9) retrieved 24 April 2024 from <u>https://medicalxpress.com/news/2021-12-head-mounted-microscope-long-term-brain-imaging.html</u>

This document is subject to copyright. Apart from any fair dealing for the purpose of private



study or research, no part may be reproduced without the written permission. The content is provided for information purposes only.