


# Rapid antigen COVID-19 tests may not be keeping pace with variants of concern, researchers find

May 5 2022, by Jacqueline Mitchell

---



Credit: Pixabay/CC0 Public Domain

By now, many have become all too familiar with rapid antigen tests for diagnosis of COVID-19. Slightly less accurate than the genetic tests

administered by healthcare professionals, the rapid antigen tests' ease-of-use allows the general public to monitor themselves for COVID-19 infections at home and make timely decisions to help stop the spread. The home tests, however, were developed and granted federal authorization for use with the original SARS-CoV-2 viral strain that emerged late in 2019. Since then, the virus has evolved countless times, and several viral variants of concern have emerged, including the highly infectious delta and omicron variants that swept across the United States last summer and winter respectively.

In a recent study, scientists at the Harvard T.H. Chan School of Public Health and Beth Israel Deaconess Medical Center used live [virus](#) culture to assess how well four rapid antigen tests are able to detect these COVID variants of concern. The findings suggest that while the rapid antigen tests remain a useful tool for the detection of COVID-19 infections, continuous assessment and updating is likely needed in the context of variants of concern. Their work is published in the *Journal of Clinical Microbiology*.

"Unlike the sensitive molecular tests that detect multiple SARS-CoV-2 genes, rapid antigen tests target a single viral protein," said co-corresponding author James Kirby, MD, director of the Clinical Microbiology Laboratory at BIDMC. "As the pandemic continues, however, some hypothesize that the performance of available antigen test may vary among the COVID variants of concern."

Using three strains of cultured live virus, the team assessed differences in the limits of detection (LoD)—the smallest amount of viral antigen detectable at 95 percent certainty—of four commercially available rapid antigen tests; the Binax, CareStart, GenBody and LumiraDx tests. The researchers found that all four tests were as sensitive to the omicron variant, if not more, as they were to original SARS-CoV-2 [viral strain](#), known as WA1. Three tests, however, showed less sensitivity to the delta

strain, with only the CareStart demonstrating equal detection of all three strains.

"We expect that the observed loss in delta sensitivity could have resulted in a 20 percent or more loss of detection in potentially infectious individuals—nevertheless, the most infectious individuals still should have been detected," said Kirby, also a professor of pathology at Harvard Medical School. "However, our findings suggest that antigen test performance needs to be reevaluated for emerging variants to ensure they still meet the intended public health testing goals of the pandemic."

**More information:** Sydney Stanley et al, Limit of Detection for Rapid Antigen Testing of the SARS-CoV-2 Omicron and Delta Variants of Concern Using Live-Virus Culture, *Journal of Clinical Microbiology* (2022). [DOI: 10.1128/jcm.00140-22](https://doi.org/10.1128/jcm.00140-22)

Provided by Beth Israel Deaconess Medical Center

Citation: Rapid antigen COVID-19 tests may not be keeping pace with variants of concern, researchers find (2022, May 5) retrieved 10 May 2024 from <https://medicalxpress.com/news/2022-05-rapid-antigen-covid-pace-variants.html>

<p>This document is subject to copyright. Apart from any fair dealing for the purpose of private study or research, no part may be reproduced without the written permission. The content is provided for information purposes only.</p>
--