


Researchers employ microfluidic diffusional sizing characterize neutralizing antibody affinity to SARS-CoV-2

March 26 2024



Schematic depict of the principles of the neutralizing antibody detection by microfluidic diffusional sizing. Credit: *Annals of Biomedical Engineering* (2024). DOI: 10.1007/s10439-024-03478-0

SARS-CoV-2 has rampantly spread around the globe and continues to cause unprecedented loss through ongoing waves of (re)infection. Increasing our understanding of the protection against infection with SARS-CoV-2 is critical to ending the pandemic.

Serological assays have been widely used to assess immune responses, but secretory antibodies, the essential first line of defense, have been studied to only a limited extent. Of particular interest and importance are neutralizing antibodies, which block the binding of the spike protein of SARS-CoV-2 to the human receptor angiotensin-converting enzyme-2 (ACE2) and thus are essential for immune defense.

In a study [published](#) in *Annals of Biomedical Engineering* researchers employed microfluidic diffusional sizing (MDS), an immobilization-free technology, to characterize neutralizing antibody affinity to SARS-CoV-2 spike receptor-binding domain (RBD) and spike trimer in [saliva](#).

Affinity measurement was obtained through a contrived sample and buffer using recombinant SARS-CoV-2 RBD and monoclonal antibody. Limited saliva samples demonstrated that MDS applies to saliva neutralizing antibody measurement. The ability to disrupt a complex of ACE2-Fc and [spike](#) trimer is shown. Using a quantitative assay on the patient sample, the researchers determined the affinity and binding site concentration of the neutralizing [antibodies](#).

More information: Cara O'Mahoney et al, Microfluidic Diffusional Sizing (MDS) Measurements of Secretory Neutralizing Antibody Affinity Against SARS-CoV-2, *Annals of Biomedical Engineering* (2024). [DOI: 10.1007/s10439-024-03478-0](https://doi.org/10.1007/s10439-024-03478-0)

Provided by Clemson University

Citation: Researchers employ microfluidic diffusional sizing characterize neutralizing antibody affinity to SARS-CoV-2 (2024, March 26) retrieved 27 April 2024 from <https://medicalxpress.com/news/2024-03-employ-microfluidic-diffusional-sizing-characterize.html>

This document is subject to copyright. Apart from any fair dealing for the purpose of private study or research, no part may be reproduced without the written permission. The content is provided for information purposes only.