


Single-nucleus RNA sequencing, droplet by droplet

28 August 2017

A condensed pedigree of single cell RNA-seqs



This family tree captures just a sampling of single-cell and single-nucleus RNA sequencing technologies that have burst onto the scene over the last eight years. Credit: Susannah M. Hamilton and Tom Ulrich, Broad Communications

Last year Broad researchers described a single-nucleus RNA sequencing method called sNuc-Seq. This system enabled researchers to study the gene expression profiles of difficult-to-isolate cell types as well as cells from archived tissues. Now a Broad-led team has overcome a key stumbling block to sNuc-Seq's widespread use: scale.

In a paper published in *Nature Methods*, postdoctoral fellows Naomi Habib, Inbal Avraham-Davidi, and Anindita Basu; core institute members Feng Zhang and Aviv Regev; and their colleagues reveal DroNc-Seq, a single-cell expression profiling technique that merges sNuc-Seq with microfluidics, allowing massively parallel measurement of gene expression in structurally-complicated tissues.

Researchers struggled in the past to study expression in neurons and other cells from complex tissues, like the brain, at the single-cell level. This was because the procedures for isolating the cells affected their RNA content and

did not always accurately capture the true proportions of the cell types present in a sample. Moreover, the procedures did not work for frozen archived tissues. sNuc-Seq bypassed those problems by using individual nuclei extracted from cells as a starting point instead.

sNuc-Seq, however, is a low-throughput technology, using 96- or 384-well plates to collect and run samples. To scale the method up to the level needed in order to efficiently study thousands of nuclei at a time, the team turned to microfluidics. Their inspiration: Drop-Seq, a single-cell RNA-seq (scRNAseq) technique that encapsulates single cells together with DNA barcoded-beads in microdroplets to greatly accelerate expression profiling experiments and reduce cost.

To test the new method's accuracy and speed, the team successfully benchmarked DroNc-Seq against Drop-Seq, sNuc-Seq, and other lower throughput scRNAseq methods using a mouse cell line and mouse brain [tissue](#). They also applied it to human brain tissue collected by the Genotype-Tissue Expression (GTEx) Project, finding that they could a) identify expression signatures unique to neurons, glial [cells](#), and other [cell types](#) in the brain (including rare types), and b) differentiate between closely related cell subtypes.

DroNc-Seq's robustness and accuracy suggest it could be a valuable addition to the stable of technologies being used as part of the Human Cell Atlas and other scRNAseq-based efforts.

More information: "Massively parallel single-nucleus RNA-seq with DroNc-seq," *Nature Methods* (2017). [DOI: 10.1038/nmeth.4407](https://doi.org/10.1038/nmeth.4407)

Provided by Broad Institute of MIT and Harvard

APA citation: Single-nucleus RNA sequencing, droplet by droplet (2017, August 28) retrieved 19 November 2019 from <https://medicalxpress.com/news/2017-08-single-nucleus-rna-sequencing-droplet.html>

This document is subject to copyright. Apart from any fair dealing for the purpose of private study or research, no part may be reproduced without the written permission. The content is provided for information purposes only.