


# Artificial intelligence may help reduce gadolinium dose in MRI

26 November 2018



Example of full-dose, 10 percent low-dose and algorithm-enhanced low-dose. Credit: Radiological Society of North America

Researchers are using artificial intelligence to reduce the dose of a contrast agent that may be left behind in the body after MRI exams, according to a study being presented today at the annual meeting of the Radiological Society of North America (RSNA).


Gadolinium is a heavy metal used in contrast material that enhances images on MRI. Recent studies have found that trace amounts of the metal remain in the bodies of people who have undergone exams with certain types of gadolinium. The effects of this deposition are not known, but radiologists are working proactively to optimize [patient safety](#) while preserving the important information that gadolinium-enhanced MRI scans provide.

"There is concrete evidence that gadolinium deposits in the brain and body," said study lead author Enhao Gong, Ph.D., researcher at Stanford University in Stanford, Calif. "While the implications of this are unclear, mitigating potential patient risks while maximizing the clinical value of the MRI exams is imperative."

Dr. Gong and colleagues at Stanford have been studying [deep learning](#) as a way to achieve this goal. Deep learning is a sophisticated artificial intelligence technique that teaches computers by examples. Through use of models called convolutional neural networks, the computer can not only recognize images but also find subtle distinctions among the imaging data that a human observer might not be capable of discerning.

To train the [deep learning algorithm](#), the researchers used MR images from 200 patients who had received contrast-enhanced MRI exams for a variety of indications. They collected three sets of images for each patient: pre-contrast scans, done prior to contrast administration and referred to as the zero-dose scans; low-dose scans, acquired after 10 percent of the standard gadolinium dose administration; and full-dose scans, acquired after 100 percent dose administration.

The algorithm learned to approximate the full-dose scans from the zero-dose and low-dose images. Neuroradiologists then evaluated the images for contrast enhancement and overall quality.



Imaging protocol used with 3 different MR series at different contrast doses. Credit: Radiological Society of North America

Results showed that the image quality was not significantly different between the low-dose,

algorithm-enhanced MR images and the full-dose, contrast-enhanced MR images. The initial results also demonstrated the potential for creating the equivalent of full-dose, contrast-enhanced MR images without any contrast agent use.

These findings suggest the method's potential for dramatically reducing gadolinium dose without sacrificing diagnostic quality, according to Dr. Gong.

"Low-dose [gadolinium images](#) yield significant untapped clinically useful information that is accessible now by using deep learning and AI," he said.

Now that the researchers have shown that the method is technically possible, they want to study it further in the clinical setting, where Dr. Gong believes it will ultimately find a home.

Future research will include evaluation of the algorithm across a broader range of MRI scanners and with different types of contrast agents.

"We're not trying to replace existing imaging technology," Dr. Gong said. "We're trying to improve it and generate more value from the existing information while looking out for the safety of our patients."

Provided by Radiological Society of North America

APA citation: Artificial intelligence may help reduce gadolinium dose in MRI (2018, November 26) retrieved 13 November 2019 from <https://medicalxpress.com/news/2018-11-artificial-intelligence-gadolinium-dose-mri.html>

*This document is subject to copyright. Apart from any fair dealing for the purpose of private study or research, no part may be reproduced without the written permission. The content is provided for information purposes only.*