

Metal stents are effective treatment for blocked bile ducts

March 20 2013

A multi-center analysis, led by Weill Cornell Medical College and published in the *Journal of Clinical Gastroenterology*, shows the use of temporary "fully covered self-expanding metal stents" (FCSEMS) can effectively fix a painful and potentially life-threatening benign biliary stricture—a severely blocked or narrowed bile duct.

"Benign biliary strictures can be managed and resolved with these new fully covered [metal stents](#) with flared ends. Our study findings are similar to the recent European study results in patients overseas who received the metal stents," says Dr. Michel Kahaleh, chief of endoscopy in the Division of Gastroenterology and Hepatology, Department of Medicine, at the Center for Advanced Digestive Care at NewYork-Presbyterian Hospital/Weill Cornell Medical Center and professor of clinical medicine at Weill Cornell Medical College.

Narrowing or blockage of a bile duct can lead to the dangerous leakage of bile into the body's organs and bloodstream. Bile is a necessary fluid produced by the liver that assists in the body's digestive process. [Bile ducts](#) are the series of thin, connected tubes that transport bile from the liver through the body to the [small intestines](#). Patients who experience a bile duct stricture may experience jaundice, severe pain, itchiness, weight-loss, elevated [liver function tests](#), cirrhosis and [liver failure](#) over a period of a few weeks.

Biliary strictures can be cancerous or benign. A benign stricture can occur as a result of surgical injury or trauma, pancreas inflammation,

scarring or hardening of bile duct tissue caused by [chronic inflammation](#), recurrent gallstones or poor reconnection of tissues post-[liver transplantation](#). A biliary stricture is often diagnosed with explorative endoscopy—a procedure that uses a thin flexible camera probe to hunt for blockages, [CT imaging](#) or MRI.

Surgery once played a pivotal role in the management of benign biliary strictures, but it's associated with significant complications and mortality. Today, minimally invasive endoscopic stent placement has emerged as a less invasive and more successful therapeutic option.

Only permanent plastic stents are currently FDA approved for benign stricture indications. However, these plastic stents involve multiple endoscopic procedures and multiple large stents. Also, plastic stents are often associated with stent blockages and failure of long-term durable results for stricture resolution.

FCSEMS are being tested in patients with benign biliary strictures to overcome the limitations of plastic stents, and the complications of their uncovered or semi-covered metal stent counterparts. The fully covered membrane design of FCSEMS prevents embedding and ingrowth of the stent into the bile duct's wall.

The research study examined 133 patient cases performed between April 2009 and October 2010 at six centers. Each patient received a FCSEMS during an endoscopic stenting procedure to open his or her severely narrowed or blocked benign bile duct stricture and restore bile flow. A 40, 60 or 80 millimeter length metal stent was implanted via a thin catheter to the blockage site over a guide wire. The stent was then partially deployed and pulled against the stricture prior to full deployment. These temporary metal stents were implanted all in one endoscopic procedure and removed after an average of three months.

This study is the first to analyze the safety and efficiency of these implanted metal stents for benign stricture resolution after their removal. "The study showed the majority of patients who received a metal stent had symptom relief, negative imaging for stricture and no need for repeat stent replacement except the patients who underwent liver transplant," says Dr. Kahaleh.

Stricture resolution after FCSEMS removal was 91.6 percent successful in post-surgical cases, 84.2 percent in gallstone related disease, chronic pancreatitis 80.7 percent, 80 percent in other etiology and, in more complex anastomotic strictures cases, it was 61.2 percent successful. Interestingly, patients who had stents placed for more than three months were 4.3 times more likely to have resolved strictures, while those with non-migrated stents were 5.4 times more likely to have resolved strictures. Also, nearly 73 percent of metal stents were successfully removed after an average of three months.

"Fully covered self-expanding metal stents are emerging as a promising therapeutic option in the treatment of benign bile duct strictures," Dr. Kahaleh says. "We found very acceptable, high rates of stricture resolution after the use and removal of metal stents in all types of patients groups. Also, those with longer duration of stent placement that didn't migrate had an increased likelihood of stricture resolution."

Current data and the FDA support the use of plastic stents in the endoscopic treatment of benign bile duct blockages. However, this current research study reports previous plastic stenting doesn't always hold up. In fact, the study shows failed plastic stenting had been performed in more than half of the study patients who underwent implant of metal stents. Also, researchers report previous plastic stenting was associated with a decreased rate of stricture resolution.

"It's time to revisit the therapeutic options we offer our patients for this

often sudden and painful benign bile duct narrowing and blockages," says Dr. Kahaleh. "Our study shows that any type of benign biliary stricture except those related to liver transplant can be treated with these metal stents, and metal stents may just do a better job fixing benign strictures than plastic stents."

The researchers report FCSEMS did have some minor complications in the patient population studied. After implantation, complications included stent migration (10.5 percent), post-procedural pain requiring hospital admission (6 percent), pancreatitis (2 percent) and stent occlusion (3 percent). Also, minimal stent removal complications included mild bleeding (1 percent), unraveling of the stent (2 percent) and hyperplastic reaction (1 percent).

According to the research team, further research investigation of FCSEMS is necessary in a prospective, multi-center randomized control trial to assess the long-term efficacy of these metal stents with a more than three-month average period of stent indwelling time.

Provided by Weill Cornell Medical College

Citation: Metal stents are effective treatment for blocked bile ducts (2013, March 20) retrieved 2 May 2024 from <https://medicalxpress.com/news/2013-03-metal-stents-effective-treatment-blocked.html>

<p>This document is subject to copyright. Apart from any fair dealing for the purpose of private study or research, no part may be reproduced without the written permission. The content is provided for information purposes only.</p>
--